

NEW

AS SUPER 4 DIGITAL

ELECTRO ACUPUNCTURE NEEDLE STIMULATOR

Features

- ❖ 30 fixed programs (including HAN)
- ❖ 7 free programs
- ❖ 4-channels – up to 8 needles
- ❖ Digital
- ❖ Sound ON/OFF
- ❖ Interlocking keypad function
- ❖ AKS for skin protection
- ❖ Intensity turns OFF if circuit is not closed for better safety
- ❖ Display indicates battery level
- ❖ Possibility to pause the stimulation
- ❖ Intensity settings in 0.1 mA steps
- ❖ Battery operation or mains supply



Technical Data

Output current: max. 20 mA

Frequency: 1-120 Hz

Pulse width: 50-300 μ s

Dimensions: 17.5 x 11.3 x 4.5 cm

Weight: 191 g (without batteries)