

TENS

EMS



	TENSTEM ECO BASIC	TENS ECO 2	SPORÉCUP XTR2	SPORÉCUP XTR4
PRE-SET PROGRAMS	12	12	30	30
FREE PROGRAMS	4	12	0	0
▀ PAIN	✓	✓	✓	✓
▀ MUSCLE STIM/REHAB	✓	✓	✓	✓
▀ VASCULAR			✓	✓
▀ SPORTS/FITNESS			✓	✓
▀ RELAXATION	✓	✓	✓	✓
▀ BEAUTY			✓	✓
CHANNELS	2	2	2	2 X 2
CARRY ON	BELT CLIP	NECK BAND, BELT CLIP	NECK BAND	NECK BAND
SOUND ON/OFF	✓	✓		
LOCKING KEY PAD	✓	✓	✓	
RECORD PATIENT DATA	✓	✓		
SKIN PROTECTION	✓	✓	✓	✓
BATTERY LEVEL	✓	✓	✓	✓
INTENSITY STEPS	1 mA	0.5 mA	0.5 mA	0.5 mA
POWER BATTERY	1.5 V AA BATTERIES	INTEGRATED BATTERY	INTEGRATED BATTERY	INTEGRATED BATTERY
POWER MAINS CIRCUIT				



	SPORÉCUP PR04	EMP 2 PRO	EMP 4 PRO	EMP 4 ECO+
PRE-SET PROGRAMS	30	21	60	25
FREE PROGRAMS	0	21	8	20
▀ PAIN	✓	✓	✓	✓
▀ MUSCLE STIM/REHAB	✓	✓	✓	✓
▀ VASCULAR	✓	✓	✓	✓
▀ SPORTS/FITNESS	✓		✓	✓
▀ RELAXATION	✓	✓	✓	✓
▀ BEAUTY	✓		✓	
CHANNELS	4	2	4	2 X 2
CARRY ON	BELT CLIP	NECK BAND	BELT CLIP	NECK BAND
SOUND ON/OFF			✓	
LOCKING KEY PAD		✓	✓	
RECORD PATIENT DATA			✓	
SKIN PROTECTION	✓	✓	✓	✓
BATTERY LEVEL	✓	✓	✓	✓
INTENSITY STEPS	1 mA	0.5 mA	1 mA	0.5 mA
POWER BATTERY	INTEGRATED BATTERY	INTEGRATED BATTERY	INTEGRATED BATTERY	INTEGRATED BATTERY
POWER MAINS CIRCUIT	✓		✓	