

TENS Overview



	TENStem eco basic	TENStem eco	TENS eco 2	EMP 2 PRO	EMP 4 PRO
PRE-SET PROGRAMS	12	12	12	21	60
PAIN	✓	✓	✓	✓	✓
MUSCLE STIM / REHAB	✓	✓	✓	✓	✓
VASCULAR				✓	✓
SPORTS / FITNESS				✓	✓
BEAUTY / RELAXATION	✓	✓	✓	✓	✓
FREE PROGRAMS	3	12	12	21	8
PAIN	✓	✓	✓	✓	✓
MUSCLE STIM / REHAB	✓	✓	✓	✓	✓
BEAUTY / RELAXATION	✓	✓	✓	✓	✓
CHANNELS	2	2	2	2	4
CARRY ON	belt clip	belt clip	neck band	belt clip	belt clip
SOUND ON/OFF	✓				✓
LOCKING KEY PAD	✓	✓	✓	✓	✓
RECORD PATIENT DATA	✓	✓	✓		✓
SKIN PROTECTION	✓	✓	✓	✓	✓
BATTERY LEVEL	✓	✓	✓	✓	✓
INTENSITY STEPS	1 mA	0.5 mA	0.5 mA	0.5 mA	1 mA
POWER BATTERY	1.5V AA batteries	9V battery	integrated rechargeable battery	1.5V AAA batteries	integrated rechargeable battery
POWER MAINS CIRCUIT				✓	✓